



Y S I Life Sciences

## YSI 2340 Data Logging Software Kit for YSI 2300 STAT Plus™ and YSI 2700 SELECT™ Analyzers

A Data Logging Software Kit is now available for the YSI 2300 STAT Plus Glucose & Lactate Analyzer and the YSI 2700 SELECT Biochemistry Analyzer. The LabVIEW-based software kit collects sample and calibration data, which can be exported to Microsoft Excel on a PC or a laptop. Up to four (4) YSI 2300 or 2700 instruments may be connected to the computer through serial or USB ports, using adapters.



### The YSI 2340 Software Kit features:

- ▲ Data recorded in Excel-compatible format on any computer with Microsoft Windows 2000, Windows XP, or Windows Vista operating system
- ▲ RS-232 cable included
- ▲ USB adapter and cable included

Data  
for Life.®

Call your local  
representative or  
the YSI Life Sciences  
technical support team

**800 659-8895**

+1 937 767-7241  
support@ysi.com

YSI Incorporated  
1700/1725 Brannum Lane  
Yellow Springs, Ohio 45387

Data for Life, 2300 STAT Plus, and  
2700 SELECT are trademarks of YSI  
Incorporated.

LabVIEW is a trademark of National  
Instruments in the United States and  
other countries.

Microsoft Excel and Windows are  
registered trademarks of Microsoft  
Corporation in the United States and  
other countries.

### Chemistries Available:

YSI 2700

Glucose

Lactate

Glutamate

Glutamine

Ethanol

Lactose

Sucrose

Galactose

Hydrogen Peroxide

Methanol

Starch

Choline

Xylose

YSI 2300

Glucose

Lactate